

TBI-FORENSIC SERVICES DIVISION PERFORMANCE VERIFICATION FORM

Instrument/Method to Verify:

Dobby (CN14383091/US1439L407)

Check the appropriate box: New Validation ☐ Modification ☐ Performance Check ☒

Laboratory: NCL ☐ KCL ☒ JCL ☐ Unit: Forensic Chemistry

Based on training and competency, the following analyst(s) are authorized by the Technical Leader to participate in this study:

Name Jessica Sosa Title SA/FS II Initials JS

Completed: *Jessica Sosa* 07/18/23

Name _____ Title _____ Initials _____

Completed:

Name _____ Title _____ Initials _____

Completed:

Verification of final review:

Unit Supervisor: Mollie Hood 7/18/2023

Section Manager: Adam Gray
Adam Gray (Jul 19, 2023 08:31 EDT)

Unit Supervisor:

Section Manager:

Unit Supervisor:

Section Manager:

Regional Administrator:

Kelvin Woodby
Kelvin Woodby (Jul 19, 2023 12:26 EDT)

IT Manager:

Matthew Buck

Regional Administrator:

Regional Administrator:

Quality Assurance Manager: Chad Johnson

By signing below, I have verified this study meets all requirements and I deem this procedure/instrument fit for the intended use:

Technical Leader:

Brett Trotter

All information must be submitted to the Quality Assurance Manager prior to the implementation of the procedure/instrument. If applicable, include training plan and competency testing schedule.

¹Validation Study for updating Dobby to Windows 10

Tennessee Bureau of Investigation Knoxville Crime Lab

7-17-23

1. Purpose of validation

The purpose of this validation study is to verify the performance of Dobby. This study will verify the methods, performance, and accuracy of identifications by comparative analysis. This is required as Dobby received firmware updates, and a new computer, requiring a limited revalidation to demonstrate that it still produces acceptable data for casework.

2. Objectives

- a. Demonstrate the accurate operation of the gas chromatograph mass spectrometer after receiving firmware updates to the instrument and the computer being replaced.
- b. Demonstrate the ability of the MassHunter and Chemstation software to function as intended software to continue to function with the updated system software.

3. Setup and Procedures²

Forensic chemistry MS identification methods were created with the same parameters as previously validated methods on this instrument before being updated.

Appropriate standards were run on each auto-sample method to demonstrate performance and consistency. Two different standards/standard mixes were run on each auto-sample method where appropriate. Other methods were verified by running the standards necessary for that method.

Various separate standard mixes were used depending on the method; xylazine, anilino-N-phenethylpiperidine (4-ANPP), heroin, para-fluorofentanyl and fentanyl. Another mix containing cocaine, hydrocodone, and pethidine was also used on methods where appropriate to compare. A final mix contained amphetamine, phentermine, and methamphetamine was used when appropriate to compare. A buprenorphine standard was also run on the appropriate methods, as well as testosterone propionate, diazepam, sildenafil, GHB, metonitazene, delta 8-THC and delta 9-THC.

The spectra obtained from these standards were compared to library spectra and the chromatographs were compared to chromatographic data obtained before the update.

4. Results and Conclusion

Analysis of the standards run on Dobby have consistently shown compliance with standard operating procedures in chapter 21. All methods have been utilized in the analysis of known drug standards. The data obtained from these standards is consistent with previously verified standards in the Knoxville forensic chemistry standard library. The instrument response was compared to older data taken from before the period of maintenance. The chromatography, peak shape, and separation all appear to be consistent to the data that the instrument was generating prior to the update.

5. Training and Competency

The instrument used to verify these methods is configured the same as the other existing instruments in the laboratory. There are no significant changes in instrumentation that would warrant an initial demonstration of proficiency for these methods.

6. Conclusions

After monitoring parameters specified in Forensic Chemistry Standard Operating Procedure found in chapter 21 and verifying the functionality of the methods with known drug standards, it has been demonstrated that Dobby produced comparable results to previous instrument response and GCMS libraries. This confirms the correct operation of the updated firmware and software.

Materials

Computer:

Device	Dobby
Processor	12th Gen Intel(R) Core(TM) i5-12500 3.00 GHz
Installed RAM	16.0 GB (15.7 GB usable)
Device ID	5D29E47F-AAC0-4B6A-A4C4-979E9DE7EB40
Product ID	00355-63303-88737-AAOEM
System type	64-bit operating system, x64-based processor
Pen and touch	No pen or touch input is available for this display
Edition	Windows 10 Pro
Version	22H2
Installed on	6/20/2023
OS build	19045.3208
Experience	Windows Feature Experience Pack 1000.19041.1000.0

Productivity tools:

Enhanced Chemstation

F.01.03.2357

Agilent MassHunter Workstation:

Version 10.0.368

GC:

Agilent 7890B

Serial: CN14383091

MS:

Agilent 5977A

Serial: US1439L407

Column:

MS: DB-1MS, Agilent 128-0112, 12 meter, 0.20 mm ID, 0.33 µm film thickness, max 360

Buprenorphine HCl

MeOH (Fisher lot# 197710 Bottle #2/4 exp. u/c)

(Sigma# SLBN0684V)

Made by LMM 9-8-22 exp. 9-8-23

Delta 8 THC

MeOH

(Cayman# 653352)

Made by CNS 6/14/24 exp. 7/14/25

Delta 9 THC

MeOH (from ampule^2)

Cerilliant/FE (#08221804)

Made by JDG and LMM 8/22/19 exp. 9/30/23

Diazepam

MeOH (Fisher lot# 211779 exp. u/c)

(Sigma# 105F0451)

Made by LMM 3/3/23 exp. 3/3/24

GHB Standard

MeOH (Fisher lot# 217364 Bottle 4/4)

(Sigma# 089H3410)

Made by CNS 6/14/23 exp. 6/14/24

GCIR Standard Mix Retention Time Mix

CHCl₃ (ACROS lot# B0542909B Bottle ¼ exp. u/c)

W/ DI H₂O cap + NH₄OH (BDH 3014-500 MLP, lot# 2014061276 exp. u/c)

Amphetamine (Sigma# 79H0603)

Methamphetamine (Sigma# 50H02161)

Phentermine HCl (Sigma# 074K1059)

Made by ZJ 2/2/23 exp. 2/1/24

Metonitazene

MeOH (Fisher lot# 211779)

(Cayman# 0608348-19)

Made by MJH 6/27/23 exp. 8/3/23

P-MOX Retention Time Mix

CHCl₃ (ACROS lot# B0542909B Bottle #1/4)

w/ DI H₂O + base cap NH₄OH (BDH lot# 3014-500MLP)

Phentermine (Sigma# 074K1059)

Methamphetamine (Sigma# 071K1580)

Xylazine (Cayman# 0471462-19)

Heroin (Cayman# 0559235-10)

Para-Fluorofentanyl (Cayman# 0541563-43)

Fentanyl (Cayman# 0549959-44)

4-ANPP (Cayman# 0540831-18)

Made by NAH 6/13/23 exp. 6/12/24

Sildenafil

MeOH (Fisher lot# 217364 Bottle 4/4)

(Cayman# 0464768-1)

Made by CNS 6/15/23 exp. 6/14/24

STD Mix Retention Time Mix

CHCl₃ (ACROS Fisher lot# B0542909B exp. u/c)

w/ DI H₂O + base cap NH₄OH (BDH lot# 2014042806 exp. u/c)

Pethidine (S.B. Penick# 413NMN-027)

Cocaine (Sigma# 69H1240)

Hydrocodone (Cayman Chem# 0524561-22-STD 1/2)

Hydrocodone (Cayman Chem# 0524561-22-STD 2/2)

Made by LMM 5/19/23 exp. 5/19/24

Testosterone Propionate

MeOH (Fisher lot# 223204 exp. u/c)

(Sigma# 034K0695)

Made by CES 5/24/23 exp. 5/24/24

**GC and GC/MS
methods**

Gas Chromatography /Mass Spectrometry methods

Auto-sample Injection	
AMINE	65 ⁰ -280 ⁰ C for 10.875 minutes
BENZOMS	210 ⁰ -280 ⁰ C for 7.0333 minutes
BUPMS	240 ⁰ -280 ⁰ C for 11.5 minutes
CFMS	150 ⁰ C for 4 minutes
DRUG1	80 ⁰ -280 ⁰ C for 14.667 minutes
GHB	80 ⁰ -280 ⁰ C for 14.667 minutes (solvent delay off)
LOWCON	80 ⁰ -280 ⁰ C for 14.667 minutes
OPIATE	200 ⁰ -280 ⁰ C for 6.6667 minutes
STEROID	70 ⁰ -280 ⁰ C for 20 minutes
THCCFMS	240 ⁰ -280 ⁰ C for 6.6667 minutes
THCMS	200 ⁰ -240 ⁰ C for 8 minutes
XLDRUGS	70 ⁰ -299 ⁰ C for 28.058 minutes

method	standard name	drugs in mixtures	RT
AMINE	GCIR STD	amphetamine	3.074
		phentermine	3.309
		methamphetamine	3.465
BENZOMS	DIAZEPAM STD	diazepam	2.643
	STD MIX	cocaine	2.031
		hydrocodone	2.677
BUPMS	BUPRENORPHINE STD	buprenorphine	7.429
	STD MIX	hydrocodone	1.550
DRUG1	GCIR STD	amphetamine	3.272
		phentermine	3.543
		methamphetamine	3.699
	METONITAZENE STD	metonitazene	13.402
	P-MOX STD	phentermine	3.537
		methamphetamine	3.692
		xylazine	7.112
		4-ANPP	8.651
		heroin	9.172
		para-fluorofentanyl	9.368
		fentanyl	9.503
		pethidine	6.394
	STD MIX	cocaine	7.902
		hydrocodone	8.621
GHB	GHB STD	GHB	1.208
LOWCON	STD MIX	pethidine	6.402
		cocaine	7.914
		hydrocodone	8.639
OPIATE	STD MIX	cocaine	2.830
		hydrocodone	3.561
STEROID	BUPRENORPHINE STD	buprenorphine	14.265
	P-MOX STD	phentermine	3.097
		methamphetamine	3.219
		xylazine	6.444
		4-ANPP	7.983
		heroin	8.504
		para-fluorofentanyl	8.700
		fentanyl	8.835
	TESTOSTERONE PROPIONATE STD	testosterone propionate	9.397
THCMS	DELTA-8-THC STD	delta-8-THC	5.580
	DELTA-9-THC STD	delta-9-THC	5.750
XLDRUGS	P-MOX STD	phentermine	3.144
		methamphetamine	3.339
		xylazine	10.432
		4-ANPP	11.931
		heroin	12.370
		para-fluorofentanyl	12.505
		fentanyl	12.594
	SILDENAFIL STD	sildenafil	20.890

METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\AMINE.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\AMINE.M
Mon Jul 17 11:27:02 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC

GC Summary
Run Time 10.875 min
Post Run Time 0 min

Oven

Temperature
Setpoint On
(Initial) 65 °C
Hold Time 1 min
Post Run 80 °C
Program
#1 Rate 30 °C/min
#1 Value 125 °C
#1 Hold Time 1 min
#2 Rate 40 °C/min
#2 Value 280 °C
#2 Hold Time 3 min

Equilibration Time 0.5 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS

Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec

Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap	
Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He	
Mode	Split
Heater	On 250 °C
Pressure	On 8.3067 psi
Total Flow	On 49 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	45 :1
Split Flow	45 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line)	
Temperature	
Setpoint	On
(Initial)	280 °C

Column	
Column #1	
Flow	On
Setpoint	1 mL/min
(Initial)	1.0812 mL/min
Post Run	
Column Information	Agilent 128-0112
DB-1ms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	65 °C
Pressure	8.3067 psi
Flow	1 mL/min
Average Velocity	58.044 cm/sec
Holdup Time	0.34457 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	65 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	

Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2

Pressure

Setpoint

(Initial)

Post Run

Off
10 psi
0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Save

Data Rate

Dual Injection Assignment

Test Plot
Off
50 Hz
Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\AMINE.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	2.50
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
2.50	40	550	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File : Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GCIR-AM.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 12:52 using AcqMethod AMINE.M
Sample Name: GCIR STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GCIR-AM.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 12:52 using AcqMethod AMINE.M
Sample Name: GCIR STD
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GCIR-AM.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 12:52 using AcqMethod AMINE.M
Sample Name: GCIR STD
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GCIRB-AM.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 12:37 using AcqMethod AMINE.M
Sample Name: GCIR STD BLANK
Misc Info : CHCl3



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\BENZOMS.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\BENZOMS.M
Mon Jul 17 11:29:14 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC
GC Summary
Run Time 7.0333 min
Post Run Time 0 min

Oven
Temperature
Setpoint On
(Initial) 210 °C
Hold Time 0.5 min
Post Run 80 °C
Program
#1 Rate 30 °C/min
#1 Value 280 °C
#1 Hold Time 4.2 min

Equilibration Time 0 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS
Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Split
Heater	On 250 °C
Pressure	On 20 psi
Total Flow	On 41.727 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	30.0001877747063 :1
Split Flow	37.478 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line)	
Temperature	
Setpoint	On
(Initial)	280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1.2493 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	210 °C
Pressure	20 psi
Flow	1.2493 mL/min
Average Velocity	68.687 cm/sec
Holdup Time	0.29117 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	210 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2

Pressure

Setpoint

(Initial)

Post Run

Off

10 psi

0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Test Plot

Save

Off

Data Rate

50 Hz

Dual Injection Assignment

Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENDOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\BENZOMS.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	1.50
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
1.50	40	500	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23DIAZE-BEN'.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 15:03 using AcqMethod BENZOMS.M
Sample Name: DIAZEPAM STD
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23DIAZEB-BEN'.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 14:54 using AcqMethod BENZOMS.M
Sample Name: DIAZEPAM STD BLANK
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23STD MIX-BEN'.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 14:27 using AcqMethod BENZOMS.M
Sample Name: STD MIX
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23STDMIX-BEN'.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 14:27 using AcqMethod BENZOMS.M
Sample Name: STD MIX
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23STDMIXB-BEN'.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 14:19 using AcqMethod BENZOMS.M
Sample Name: STD MIX BLANK
Misc Info : CHCl3



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\BUPMS.M

Method Sections To Run:

() Save Copy of Method With Data
() Instrument Control Pre-Run Cmd/Macro =
() Data Analysis Pre-Run Cmd/Macro =
(X) Data Acquisition
(X) Data Analysis
() Instrument Control Post-Run Cmd/Macro =
() Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\BUPMS.M

Mon Jul 17 11:30:01 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC

GC Summary

Run Time 11.5 min
Post Run Time 0 min

Oven

Temperature
Setpoint On
(Initial) 240 °C
Hold Time 0 min
Post Run 80 °C
Program
#1 Rate 40 °C/min
#1 Value 280 °C
#1 Hold Time 10.5 min

Equilibration Time 0.5 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS

Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Split
Heater	On 250 °C
Pressure	On 17.958 psi
Total Flow	On 74 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	69.9998743249809 :1
Split Flow	70 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line) Temperature	
Setpoint (Initial)	On 280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	240 °C
Pressure	17.958 psi
Flow	1 mL/min
Average Velocity	62.05 cm/sec
Holdup Time	0.32232 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	240 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2
Pressure
Setpoint
(Initial)
Post Run

Off
10 psi
0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Test Plot

Save

Off

Data Rate

50 Hz

Dual Injection Assignment

Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\BUPMS.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:30:06 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\BUPMS.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	1.00
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
1.00	40	500	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23BUP-BP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 16:50 using AcqMethod BUPMS.M
Sample Name: BUPRENORPHINE STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23BUPB-BP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 16:36 using AcqMethod BUPMS.M
Sample Name: BUPRENORPHINE STD BLANK
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIX-BP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 16:03 using AcqMethod BUPMS.M
Sample Name: STD MIX
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIXB-BP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 15:50 using AcqMethod BUPMS.M
Sample Name: STD MIX BLANK
Misc Info : CHCl3



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\DRUG1.M

Method Sections To Run:

() Save Copy of Method With Data
() Instrument Control Pre-Run Cmd/Macro =
() Data Analysis Pre-Run Cmd/Macro =
(X) Data Acquisition
(X) Data Analysis
() Instrument Control Post-Run Cmd/Macro =
() Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\DRUG1.M
Mon Jul 17 11:32:55 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC
GC Summary
Run Time 14.667 min
Post Run Time 0 min

Oven
Temperature
Setpoint On
(Initial) 80 °C
Hold Time 2 min
Post Run 50 °C
Program
#1 Rate 30 °C/min
#1 Value 280 °C
#1 Hold Time 6 min

Equilibration Time 0.5 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS
Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Split
Heater	On 250 °C
Pressure	On 9.1608 psi
Total Flow	On 49 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	45 :1
Split Flow	45 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line) Temperature	
Setpoint	On
(Initial)	280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	80 °C
Pressure	9.1608 psi
Flow	1 mL/min
Average Velocity	58.448 cm/sec
Holdup Time	0.34218 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	80 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2
Pressure
Setpoint
(Initial)
Post Run

Off
10 psi
0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Test Plot

Save

Off

Data Rate

50 Hz

Dual Injection Assignment

Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\DRUG1.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:33:00 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\DRUG1.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	2.50
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
2.50	40	550	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-11-23\JS-07-11-23-GCIR-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 11 Jul 2023 18:28 using AcqMethod DRUG1.M
Sample Name: GCIR STD MIX
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-11-23\JS-07-11-23-GCIR-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 11 Jul 2023 18:28 using AcqMethod DRUG1.M
Sample Name: GCIR STD MIX
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-11-23\JS-07-11-23-GCIR-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 11 Jul 2023 18:28 using AcqMethod DRUG1.M
Sample Name: GCIR STD MIX
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-11-23\JS-07-11-23-GCIR-DB.D
Operator : Dobby
Instrument : Dobby
Acquired : 11 Jul 2023 18:10 using AcqMethod DRUG1.M
Sample Name: GCIR STD MIX BLANK
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23METO-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 14:16 using AcqMethod DRUG1.M
Sample Name: METONITAZENE STD
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23METO-B.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:57 using AcqMethod DRUG1.M
Sample Name: METONITAZENE STD BLANK
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3




File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3




File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:21 using AcqMethod DRUG1.M
Sample Name: P-MOX STD
Misc Info : CHCl3




File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOXB-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 13:03 using AcqMethod DRUG1.M
Sample Name: P-MOX STD BLANK
Misc Info : CHCl3





Scan 540 (5.735 min): JS-07-13-23P-MOXB-D.D\data.ms





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 12:27 using AcqMethod DRUG1.M
Sample Name: STD MIX
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 12:27 using AcqMethod DRUG1.M
Sample Name: STD MIX
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIX-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 12:27 using AcqMethod DRUG1.M
Sample Name: STD MIX
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIXB-D.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 12:09 using AcqMethod DRUG1.M
Sample Name: STD MIX BLANK
Misc Info : CHCl3



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\GHB.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\GHB.M

Mon Jul 17 11:34:53 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC
GC Summary
Run Time 14.667 min
Post Run Time 0 min

Oven
Temperature
Setpoint On
(Initial) 80 °C
Hold Time 1 min
Post Run 80 °C
Program
#1 Rate 30 °C/min
#1 Value 280 °C
#1 Hold Time 7 min

Equilibration Time 0 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS
Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Split
Heater	On 250 °C
Pressure	On 9.1608 psi
Total Flow	On 44 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	39.9999274271017 :1
Split Flow	40 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line)	
Temperature	
Setpoint	On
(Initial)	280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	80 °C
Pressure	9.1608 psi
Flow	1 mL/min
Average Velocity	58.448 cm/sec
Holdup Time	0.34218 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	80 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2

Pressure

Setpoint

(Initial)

Post Run

Off
10 psi
0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Save

Data Rate

Dual Injection Assignment

Test Plot
Off
50 Hz
Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENDOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\GHB.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:34:59 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\GHB.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	0.00
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
0.00	40	500	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GHB.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 13:41 using AcqMethod GHB.M
Sample Name: GHB STD
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GHB.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 13:41 using AcqMethod GHB.M
Sample Name: GHB STD
Misc Info : MeOH



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GHB-B.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 13:23 using AcqMethod GHB.M
Sample Name: GHB STD BLANK
Misc Info : MeOH



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23GHB-B.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 13:23 using AcqMethod GHB.M
Sample Name: GHB STD BLANK
Misc Info : MeOH



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\LOWCON.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\LOWCON.M

Mon Jul 17 11:36:08 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC
GC Summary
Run Time 14.667 min
Post Run Time 0 min

Oven
Temperature
Setpoint On
(Initial) 80 °C
Hold Time 2 min
Post Run 50 °C
Program
#1 Rate 30 °C/min
#1 Value 280 °C
#1 Hold Time 6 min

Equilibration Time 0.5 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS
Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Splitless
Heater	On 250 °C
Pressure	On 9.1608 psi
Total Flow	On 46.7 mL/min
Septum Purge Flow	On 3 mL/min
Septum Purge Flow Mode	Standard
Gas Saver	On 20 After 2 min mL/min
Purge Flow to Split Vent	42.7 mL/min at 1 min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line) Temperature	
Setpoint	On
(Initial)	280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-1ms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	80 °C
Pressure	9.1608 psi
Flow	1 mL/min
Average Velocity	58.448 cm/sec
Holdup Time	0.34218 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	80 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2
Pressure
Setpoint
(Initial)
Post Run

Off
10 psi
0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Test Plot

Save

Off

Data Rate

50 Hz

Dual Injection Assignment

Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\LOWCON.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:36:13 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\LOWCON.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	2.50
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
2.50	40	425	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23STD MIXDIL-LOW.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 09:19 using AcqMethod LOWCON.M
Sample Name: STD MIX DILUTED
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23STDMIXDIL-LOW.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 09:19 using AcqMethod LOWCON.M
Sample Name: STD MIX DILUTED
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23STD MIXDIL-LOW.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 09:19 using AcqMethod LOWCON.M
Sample Name: STD MIX DILUTED
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23STD MIXDILB-LOW.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 09:01 using AcqMethod LOWCON.M
Sample Name: STD MIX DILUTED BLANK
Misc Info : CHCl3



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\OPIATE.M

Method Sections To Run:

() Save Copy of Method With Data
() Instrument Control Pre-Run Cmd/Macro =
() Data Analysis Pre-Run Cmd/Macro =
(X) Data Acquisition
(X) Data Analysis
() Instrument Control Post-Run Cmd/Macro =
() Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\OPIATE.M

Mon Jul 17 11:37:03 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC

GC Summary

Run Time 6.6667 min
Post Run Time 0 min

Oven

Temperature

Setpoint On
(Initial) 200 °C
Hold Time 1 min
Post Run 50 °C
Program
#1 Rate 30 °C/min
#1 Value 280 °C
#1 Hold Time 3 min

Equilibration Time 0 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS

Front Injector

Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B

Sample Overlap Mode	Sample overlap is not enabled
---------------------	-------------------------------

ALS Errors	Pause for user interaction
------------	----------------------------

Front SS Inlet He Mode	Split
Heater	On 250 °C
Pressure	On 15.806 psi
Total Flow	On 94 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	90.00000001 :1
Split Flow	90 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv

Back SS Inlet He Excluded from Readiness Mode	***Excluded from Affecting GC's Readiness State***
Heater	Split
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	Off
Split Flow	1 :1
Liner	0 mL/min
	A Liner has not been selected.

Thermal Aux 2 (MSD Transfer Line) Temperature Setpoint (Initial)	On 280 °C
--	-----------

Column
Column #1
Flow

Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	200 °C
Pressure	15.806 psi
Flow	1 mL/min
Average Velocity	61.249 cm/sec
Holdup Time	0.32654 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	200 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2

Pressure

Setpoint

(Initial)

Post Run

Off

10 psi

0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Test Plot

Save

Off

Data Rate

50 Hz

Dual Injection Assignment

Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\OPIATE.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No

Printer: Yes

File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No

Printer: Yes

File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:37:08 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\OPIATE.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	2.50
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
2.50	40	500	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIX-OP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 15:22 using AcqMethod OPIATE.M
Sample Name: STD MIX
Misc Info : CHC13



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STD MIX-OP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 15:22 using AcqMethod OPIATE.M
Sample Name: STD MIX
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23STDMIXB-OP.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 15:13 using AcqMethod OPIATE.M
Sample Name: STD MIX BLANK
Misc Info : CHCl3



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\STEROID.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\STEROID.M

Mon Jul 17 11:38:01 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC

GC Summary

Run Time 20 min
Post Run Time 0 min

Oven

Temperature

Setpoint On
(Initial) 70 °C
Hold Time 1 min
Post Run 50 °C
Program
#1 Rate 30 °C/min
#1 Value 280 °C
#1 Hold Time 12 min

Equilibration Time 0.5 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS

Front Injector

Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Split
Heater	On 250 °C
Pressure	On 8.592 psi
Total Flow	On 46 mL/min
Septum Purge Flow	Off
Gas Saver	On 20 After 2 min mL/min
Split Ratio	45 :1
Split Flow	45 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line) Temperature	
Setpoint	On
(Initial)	280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	70 °C
Pressure	8.592 psi
Flow	1 mL/min
Average Velocity	58.18 cm/sec
Holdup Time	0.34376 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	70 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2

Pressure

Setpoint

(Initial)

Post Run

Off

10 psi

0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot

Description

Test Plot

Save

Off

Data Rate

50 Hz

Dual Injection Assignment

Front Sample

Signal #2:

Description

None

Signal #3:

Description

None

Signal #4:

Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENDOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\STEROID.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No

Printer: Yes

File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No

Printer: Yes

File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:38:07 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\STEROID.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	2.50
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
2.50	40	550	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23BUP-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 18:53 using AcqMethod STER0ID.M
Sample Name: BUPRENORPHINE STD
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23BUPB-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 18:29 using AcqMethod STERIOD.M
Sample Name: BUPRENORPHINE STD BLANK
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STERIOD.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : •Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STERIOD.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STER0ID.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STEROID.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STERIOD.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STEROID.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STEROID.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STEROID.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOX-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:47 using AcqMethod STEROID.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23P-MOXB-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 17:24 using AcqMethod STEROID.M
Sample Name: P-MOX STD BLANK
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23TESTPROP-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 19:58 using AcqMethod STEROID.M
Sample Name: TESTOSTERONE PROPIONATE STD
Misc Info : MeOH



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-13-23\JS-07-13-23TESTPROPB-S.D
Operator : Dobby
Instrument : Dobby
Acquired : 13 Jul 2023 19:35 using AcqMethod STER0ID.M
Sample Name: TESTOSTERONE PROPIONATE STD BLANK
Misc Info : MeOH



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\THCMS.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\THCMS.M

Mon Jul 17 11:39:28 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC

GC Summary

Run Time 8 min
Post Run Time 0 min

Oven

Temperature

Setpoint On
(Initial) 200 °C
Hold Time 2 min
Post Run 50 °C
Program
#1 Rate 10 °C/min
#1 Value 240 °C
#1 Hold Time 2 min

Equilibration Time 0.5 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS

Front Injector

Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 2
Solvent A Washes (PostInj) 2
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 2
Solvent B Washes (PostInj) 2
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 2
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec
Sample Depth Disabled
Injection Type Standard
L1 Airgap 0.2 µL

Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He Mode	Split
Heater	On 280 °C
Pressure	On 26.206 psi
Total Flow	On 83.918 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	44 :1
Split Flow	79.12 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He Mode	***Excluded from Affecting GC's Readiness State***
Excluded from Readiness Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line)	
Temperature	
Setpoint	On
(Initial)	280 °C
Column	
Column #1	
Flow	

Setpoint	On
(Initial)	1.7982 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	200 °C
Pressure	26.206 psi
Flow	1.7982 mL/min
Average Velocity	82.132 cm/sec
Holdup Time	0.24351 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	200 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Aux EPC 3 N2
Pressure
Setpoint
(Initial)
Post Run

Off
10 psi
0 psi

Excluded from Readiness

Excluded from Affecting GC's Readiness State

Signals

Signal #1: Test Plot
Description
Save
Data Rate
Dual Injection Assignment

Test Plot
Off
50 Hz
Front Sample

Signal #2:
Description

None

Signal #3:
Description

None

Signal #4:
Description

None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\THCMS.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:39:33 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\THCMS.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	1.00
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
1.00	40	550	70	3,125 [N=1]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23DELTA-8-THCMS.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 10:46 using AcqMethod THCMS.M
Sample Name: DELTA-8-THC STD
Misc Info : MeOH





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23DELTA-8B-THCMS.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 10:35 using AcqMethod THCMS.M
Sample Name: DELTA-8-THC STD BLANK
Misc Info : MeOH



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23DELTA-9-THCMS.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 10:15 using AcqMethod THCMS.M
Sample Name: DELTA-9-THC ~~BLANK~~ *STD*
Misc Info : MeOH



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-17-23\JS-07-17-23DELTA-9B-THCMS.D
Operator : Dobby
Instrument : Dobby
Acquired : 17 Jul 2023 10:05 using AcqMethod THCMS.M
Sample Name: DELTA-9-THC STD BLANK
Misc Info : MeOH



METHOD CONTROL PARAMETERS

Method Information for: C:\MassHunter\GCMS\1\methods\XLDRUGS.M

Method Sections To Run:

- () Save Copy of Method With Data
- () Instrument Control Pre-Run Cmd/Macro =
- () Data Analysis Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Instrument Control Post-Run Cmd/Macro =
- () Data Analysis Post-Run Cmd/Macro =

Method Comments:

This is the default method

END OF METHOD CONTROL PARAMETERS

C:\MassHunter\GCMS\1\methods\XLDRUGS.M

Mon Jul 17 11:40:34 2023

Control Information

Sample Inlet : GC
Injection Source : GC ALS
Injection Location: Front
Mass Spectrometer : Enabled

No Sample Prep method has been assigned to this method.

GC

GC Summary

Run Time 28.058 min
Post Run Time 0 min

Oven

Temperature
Setpoint On
(Initial) 70 °C
Hold Time 0 min
Post Run 80 °C
Program
#1 Rate 15 °C/min
#1 Value 150 °C
#1 Hold Time 3 min
#2 Rate 40 °C/min
#2 Value 299 °C
#2 Hold Time 16 min

Equilibration Time 0 min
Max Temperature 325 °C
Maximum Temperature Override Disabled
Slow Fan Disabled

ALS

Front Injector
Syringe Size 10 µL
Injection Volume 1 µL
Solvent A Washes (PreInj) 1
Solvent A Washes (PostInj) 1
Solvent A Volume 8 µL
Solvent B Washes (PreInj) 1
Solvent B Washes (PostInj) 1
Solvent B Volume 8 µL
Sample Washes 1
Sample Wash Volume 8 µL
Sample Pumps 1
Dwell Time (PreInj) 0 min
Dwell Time (PostInj) 0 min
Solvent Wash Draw Speed 300 µL/min
Solvent Wash Dispense Speed 6000 µL/min
Sample Wash Draw Speed 300 µL/min
Sample Wash Dispense Speed 6000 µL/min
Injection Dispense Speed 6000 µL/min
Viscosity Delay 0 sec

Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Back Injector	
Syringe Size	10 µL
Injection Volume	1 µL
Solvent A Washes (PreInj)	1
Solvent A Washes (PostInj)	1
Solvent A Volume	8 µL
Solvent B Washes (PreInj)	1
Solvent B Washes (PostInj)	1
Solvent B Volume	8 µL
Sample Washes	1
Sample Wash Volume	8 µL
Sample Pumps	1
Dwell Time (PreInj)	0 min
Dwell Time (PostInj)	0 min
Solvent Wash Draw Speed	300 µL/min
Solvent Wash Dispense Speed	6000 µL/min
Sample Wash Draw Speed	300 µL/min
Sample Wash Dispense Speed	6000 µL/min
Injection Dispense Speed	6000 µL/min
Viscosity Delay	0 sec
Sample Depth	Disabled
Injection Type	Standard
L1 Airgap	0.2 µL
Solvent Wash Mode	A, B
Sample Overlap	
Mode	Sample overlap is not enabled
ALS Errors	Pause for user interaction
Front SS Inlet He	
Mode	Split
Heater	On 250 °C
Pressure	On 8.592 psi
Total Flow	On 44 mL/min
Septum Purge Flow	On 3 mL/min
Gas Saver	On 20 After 2 min mL/min
Split Ratio	39.9999274271017 :1
Split Flow	40 mL/min
Liner	Agilent 5190-2270: 900 µL (Splitless, single taper, deactiv
Back SS Inlet He	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Mode	Split
Heater	Off
Pressure	Off
Total Flow	Off
Septum Purge Flow	Off
Gas Saver	Off
Split Ratio	1 :1
Split Flow	0 mL/min
Liner	A Liner has not been selected.
Thermal Aux 2 (MSD Transfer Line)	
Temperature	
Setpoint	On
(Initial)	280 °C

Column	
Column #1	
Flow	
Setpoint	On
(Initial)	1 mL/min
Post Run	1.0812 mL/min
Column Information	Agilent 128-0112
DB-lms	
Temperature Range	-60 °C-340 °C (350 °C)
Dimensions	12 m x 200 µm x 0.33 µm
Column lock	Unlocked
In	Front SS Inlet He
Out	MSD
(Initial)	70 °C
Pressure	8.592 psi
Flow	1 mL/min
Average Velocity	58.18 cm/sec
Holdup Time	0.34376 min
Control Mode	Constant Flow
Column #2	
Pressure	
Setpoint	Off
(Initial)	0 psi
Post Run	0 psi
Column Information	Agilent 19091Z-211
HP-1	
Temperature Range	-60 °C-325 °C (350 °C)
Dimensions	15 m x 320 µm x 1 µm
Column lock	Unlocked
In	Back SS Inlet He
Out	Front Detector FID
(Initial)	70 °C
Pressure	0 psi
Flow	0 mL/min
Average Velocity	0 cm/sec
Holdup Time	∞ min
Control Mode	Constant Pressure
Column Outlet Pressure	0 psi
Front Detector FID	
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Makeup	N2
Heater	Off
H2 Flow	Off
Air Flow	Off
Makeup Flow	Off
Carrier Gas Flow Correction	Constant Makeup and Fuel Flow
Flame	Off
Aux EPC 1,2,3	
Aux EPC 1 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi
Excluded from Readiness	***Excluded from Affecting GC's Readiness State***
Aux EPC 2 N2	
Pressure	

Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness ***Excluded from Affecting GC's Readiness State***

Aux EPC 3 N2	
Pressure	
Setpoint	Off
(Initial)	10 psi
Post Run	0 psi

Excluded from Readiness ***Excluded from Affecting GC's Readiness State***

Signals	
Signal #1: Test Plot	
Description	Test Plot
Save	Off
Data Rate	50 Hz
Dual Injection Assignment	Front Sample

Signal #2:	
Description	None

Signal #3:	
Description	None

Signal #4:	
Description	None

TUNE PARAMETERS for SN: US1439L407

Trace Ion Detection is OFF.

34.593	:	EMISSION	
70.007	:	ENERGY	
34.899	:	REPELLER	
90.331	:	IONFOCUS	
20.181	:	ENTRANCE_LENS	
963.560	:	EMVOLTS	
		1151.3	: Actual EMV
		2.60	: GAIN FACTOR
1322.000	:	AMUGAIN	
123.750	:	AMUOFFSET	
1.000	:	FILAMENT	
0.000	:	DCPOLARITY	
15.312	:	ENTLENSOFFSET	
0.000	:	Ion_Body	
0.000	:	EXTLENS	
-468.000	:	MASSGAIN	
-36.000	:	MASSOFFSET	

END OF TUNE PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MassHunter\GCMS\1\methods\XLDRUGS.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No

Printer: Yes

File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search	Minimum Quality
DEMO.L	0

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No

Printer: Yes

File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Calibration Last Updated:

Reference Window: 10.00 Percent
Non-Reference Window: 5.00 Percent
Correlation Window: 0.02 minutes
Default Multiplier: 1.00
Default Sample Concentration: 0.00

Compound Information

*** Empty Quantitation Database ***

END OF DATA ANALYSIS PARAMETERS

Mon Jul 17 11:40:40 2023

Single Quadrupole Acquisition Method - MS Parameters Report

Method file	C:\MassHunter\GCMS\1\methods\XLDRUGS.M
Tune file	ATUNE.U
Ion source	EI
Source temperature (°C)	230
Quad temperature (°C)	150
Fixed Electron energy (eV)	70.0
Acquisition Type	Scan
Stop time (min)	10.00
Solvent delay (min)	2.00
Trace Ion Detection	False
Gain Factor	2.6
EM Saver	False
EM Saver Limit	N/A

Scan Time Segments

Time	Start Mass	End Mass	Threshold	Scan Speed
2.00	40	500	150	1,562 [N=2]

Timed Events

Time	Type of Event	Parameter
------	---------------	-----------



Real-Time Plots

Type of Plot	Label	Low Mass	High Mass
Total Ion	N/A	N/A	N/A
Spectrum	N/A	N/A	N/A
Extracted Ion	Scan 1-1	50	550
Extracted Ion	Scan 1-2	40	550



Self-Cleaning Ion Source Parameters

Mode	No Cleaning
------	-------------



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOX-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 10:26 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD
Misc Info : CHCl3





File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23P-MOXB-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 09:55 using AcqMethod XLDRUGS.M
Sample Name: P-MOX STD BLANK
Misc Info : CHCl3



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
... ck\Test\7-14-23\JS-07-14-23SILD-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 11:47 using AcqMethod XLDRUGS.M
Sample Name: SILDENAFIL STD
Misc Info : MeOH



File :Z:\Forensic Chemistry\Instruments\GCMS\Dobby\Performance Che
ck\Test\7-14-23\JS-07-14-23SILDB-XL.D
Operator : Dobby
Instrument : Dobby
Acquired : 14 Jul 2023 11:16 using AcqMethod XLDRUGS.M
Sample Name: SILDENAFIL STD BLANK
Misc Info : MeOH




Autotune - 5977

Tune timestamp: 6/21/2023 11:18 AM (UTC-07:00)

Dobby

C:\MASSHUNTER\GCMS\1\5977\atune.u

US1439L407




Ion Polarity	Pos	PFTBA	Open
Emission	34.6	Mass Gain	-469
Electron Energy	70.0	Mass Offset	-33
Filament	1	Amu Gain	1323
Repeller	34.90	Amu Offset	122.75
Ion Focus	87.4	Width219	-0.006
Entrance Lens	25.2	DC Polarity	Pos
Ent Lens Offset	12.90	HED Enable	On
Ion Body	0.00	EM Volts	986.4
Post Extractor 1	0	Extractor Lens	0.00
Post Extractor 2	0	Scan Speed	3
JetClean Flow Actual/[Setpoint]	0.00 [0.00]	Averages	3

Actual m/z	Abund	Rel Abund	Pw50
69.00	452,508	100.0%	0.60
218.90	211,805	46.8%	0.60
502.00	6,472	1.4%	0.62

Temperatures and Pressures		
MS Source	230 Turbo Speed	100.0
MS Quad	150 Hi Vac	8.57e-06

Low	High	Step	Speed	Threshold	Peaks	Base	Abundance	Total Ion
10.00	701.00	0.10	3	100	104	69.00	434,176	1,167,785



Target m/z	Actual m/z	Abund	Rel Abund	Iso m/z	Iso Abund	Iso Ratio
69.00	69.00	434,176	100.0%	70.00	4,479	1.0%
219.00	219.00	201,344	46.4%	220.00	8,493	4.2%
502.00	502.00	6,105	1.4%	503.00	599	9.8%

Air/Water Check: H2O ~0.4% N2 ~0.6% O2 ~0.2% CO2 ~0.1% N2/H2O ~142.3%

Column(1) Flow: 0.68 Column(2): 0.08 ml/min Interface Temp: 280

Ramp Criteria:

Ion Focus maximum 90 volts using ion 502; Electron Multiplier Gain 40460.083

Repeller maximum 35 volts using ion 219; Gain Factor 0.4046

Mass Gain Values(Scan Speed): -463(3) -456(2) -440(1) -404(0) -317(FS1) -326(FS2)

TARGET MASS:	50	69	131	219	414	502	1050
Amu Offset	122.8	122.8	122.8	122.8	122.8	122.8	122.8
Entrance Lens Offset	12.9	12.9	12.9	12.9	12.9	12.9	12.9


Autotune - 5977

Tune timestamp: 7/6/2023 12:52 PM (UTC-07:00)

C:\MASSHUNTER\GCMS\1\5977\ATUNE.U

Dobby

US1439L407




Ion Polarity	Pos	PFTBA	Open
Emission	34.6	Mass Gain	-459
Electron Energy	70.0	Mass Offset	-35
Filament	1	Amu Gain	1320
Repeller	34.90	Amu Offset	124.50
Ion Focus	90.3	Width219	0.001
Entrance Lens	20.2	DC Polarity	Pos
Ent Lens Offset	15.31	HED Enable	On
Ion Body	0.00	EM Volts	959.5
Post Extractor 1	0	Extractor Lens	0.00
Post Extractor 2	0	Scan Speed	3
JetClean Flow Actual/[Setpoint]	0.00 [0.00]	Averages	3

Actual m/z	Abund	Rel Abund	Pw50
69.00	528,182	100.0%	0.60
219.00	347,085	65.7%	0.60
502.10	16,460	3.1%	0.60

Temperatures and Pressures			
MS Source	230	Turbo Speed	100.0
MS Quad	150	Hi Vac	1.08e-05

Low	High	Step	Speed	Threshold	Peaks	Base	Abundance	Total Ion
10.00	701.00	0.10	3	100	112	69.00	502,912	1,580,144



Target m/z	Actual m/z	Abund	Rel Abund	Iso m/z	Iso Abund	Iso Ratio
69.00	69.00	502,912	100.0%	70.10	5,295	1.1%
219.00	219.00	326,848	65.0%	220.00	14,125	4.3%
502.00	502.00	15,703	3.1%	503.10	1,553	9.9%

Air/Water Check: H2O ~0.3% N2 ~0.3% O2 ~0.1% CO2 ~0.1% N2/H2O ~116.1%

Column(1) Flow: 1.00 Column(2): 0.08 ml/min Interface Temp: 280

Ramp Criteria:

Ion Focus maximum 90 volts using ion 502; Electron Multiplier Gain 28195.306

Repeller maximum 35 volts using ion 219; Gain Factor 0.2820

Mass Gain Values(Scan Speed): -455(3) -448(2) -432(1) -402(0) -324(FS1) -323(FS2)

TARGET MASS:	50	69	131	219	414	502	1050
Amu Offset	124.5	124.5	124.5	124.5	124.5	124.5	124.5
Entrance Lens Offset	15.3	15.3	15.3	15.3	15.3	15.3	15.3


Autotune - 5977

Tune timestamp: 7/10/2023 4:03 PM (UTC-04:00)

C:\MASSHUNTER\GCMS\1\5977\ATUNE.U

Dobby

US1439L407




Ion Polarity	Pos	PFTBA	Open
Emission	34.6	Mass Gain	-459
Electron Energy	70.0	Mass Offset	-36
Filament	1	Amu Gain	1319
Repeller	34.90	Amu Offset	123.94
Ion Focus	90.3	Width219	-0.007
Entrance Lens	22.7	DC Polarity	Pos
Ent Lens Offset	15.15	HED Enable	On
Ion Body	0.00	EM Volts	948.2
Post Extractor 1	0	Extractor Lens	0.00
Post Extractor 2	0	Scan Speed	3
JetClean Flow	0.00	Averages	3
Actual/[Setpoint]	[0.00]		

Actual m/z	Abund	Rel Abund	Pw50
69.00	515,373	100.0%	0.60
219.00	326,869	63.4%	0.60
502.00	15,320	3.0%	0.60

Temperatures and Pressures			
MS Source	230	Turbo Speed	100.0
MS Quad	150	HI Vac	1.12e-05

Low	High	Step	Speed	Threshold	Peaks	Base	Abundance	Total Ion
10.00	701.00	0.10	3	100	109	69.00	490,752	1,523,337



Target m/z	Actual m/z	Abund	Rel Abund	Iso m/z	Iso Abund	Iso Ratio
69.00	69.00	490,752	100.0%	70.00	5,360	1.1%
219.00	219.00	312,832	63.7%	220.00	13,428	4.3%
502.00	502.00	14,267	2.9%	503.00	1,400	9.8%

Air/Water Check: H2O ~0.3% N2 ~0.4% O2 ~0.2% CO2 ~0.1% N2/H2O ~146.4%

Column(1) Flow: 1.00 Column(2): 0.07 ml/min Interface Temp: 280

Ramp Criteria:

Ion Focus maximum 90 volts using Ion 502; Electron Multiplier Gain 24306.271

Repeller maximum 35 volts using Ion 219; Gain Factor 0.2431

Mass Gain Values(Scan Speed): -455(3) -448(2) -432(1) -402(0) -324(FS1) -323(FS2)

TARGET MASS:	50	69	131	219	414	502	1050
Amu Offset	123.9	123.9	123.9	123.9	123.9	123.9	123.9
Entrance Lens Offset	15.1	15.1	15.1	15.1	15.1	15.1	15.1


Autotune - 5977

Tune timestamp: 7/17/2023 7:16 AM (UTC-04:00)

C:\MASSHUNTER\GCMS\1\5977\ATUNE.U

Dobby

US1439L407




Ion Polarity	Pos	PFTBA	Open
Emission	34.6	Mass Gain	-468
Electron Energy	70.0	Mass Offset	-36
Filament	1	Amu Gain	1322
Repeller	34.90	Amu Offset	123.75
Ion Focus	90.3	Width219	-0.001
Entrance Lens	20.2	DC Polarity	Pos
Ent Lens Offset	15.31	HED Enable	On
Ion Body	0.00	EM Volts	963.6
Post Extractor 1	0	Extractor Lens	0.00
Post Extractor 2	0	Scan Speed	3
JetClean Flow Actual/[Setpoint]	0.00 [0.00]	Averages	3

Actual m/z	Abund	Rel Abund	Pw50
69.00	545,071	100.0%	0.60
218.90	353,205	64.8%	0.60
502.00	17,103	3.1%	0.60

Temperatures and Pressures			
MS Source	230	Turbo Speed	100.0
MS Quad	150	Hi Vac	1.07e-05

Low	High	Step	Speed	Threshold	Peaks	Base	Abundance	Total Ion
10.00	701.00	0.10	3	100	119	69.00	522,176	1,627,037



Target m/z	Actual m/z	Abund	Rel Abund	Iso m/z	Iso Abund	Iso Ratio
69.00	69.00	522,176	100.0%	70.00	5,781	1.1%
219.00	219.00	333,312	63.8%	220.00	14,383	4.3%
502.00	502.00	15,973	3.1%	503.00	1,624	10.2%

Air/Water Check: H2O ~0.3% N2 ~0.9% O2 ~0.4% CO2 ~0.1% N2/H2O ~305.0%

Column(1) Flow: 1.00 Column(2): 0.07 ml/min Interface Temp: 280

Ramp Criteria:

Ion Focus maximum 90 volts using ion 502; Electron Multiplier Gain 28341.880

Repeller maximum 35 volts using ion 219; Gain Factor 0.2834

Mass Gain Values(Scan Speed): -457(3) -450(2) -434(1) -403(0) -315(FS1) -367(FS2)

TARGET MASS:	50	69	131	219	414	502	1050
Amu Offset	123.8	123.8	123.8	123.8	123.8	123.8	123.8
Entrance Lens Offset	15.3	15.3	15.3	15.3	15.3	15.3	15.3