



TENNESSEE BUREAU OF INVESTIGATION

Forensic Services Division

Barometer Maintenance Procedure

Barometers utilized within the Breath Alcohol Calibration Unit require the following maintenance:

- The devices will be routinely checked to ensure they are properly functioning. If found inoperable, they will be taken out of service;
- The devices will be calibrated by an outside vendor every two years. Specifications for the calibration will be communicated to the vendor. An example of the vendor spec letter is available in the appendix;
- Following calibration, the barometer will be checked for consistency with another device currently in service;
- A list of barometers will be maintained to include serial number and calibration schedule.