

TENNESSEE BUREAU OF INVESTIGATION
Forensic Services Division

Forensic Chemistry Standard Operating Procedure Manual
Refrigerators, Freezers, and Ovens



6.0 REFRIGERATORS, FREEZERS, AND OVENS

Refrigerators, freezers, and ovens will be checked on a weekly basis with a NIST-traceable thermometer to ensure proper function. The temperatures will be recorded in the appropriate logbook. Log sheets are located in the Lab Documentation folder in Ensur.

NIST-traceable thermometers shall be purchased or recertified annually.

Note: Ovens that are only used to dry down samples do not need to be checked.